

User Manual

REALQUALITY

RS-MYCO P

v.2.1

code RQ-S47

Kit for detection
and quantification of
Mycoplasma pneumoniae (MYCO P)
by Real-Time PCR



1	PRODUCT INFORMATION	3	11.3 Instrument programming	18
1.1	Intended use	3	11.3.1 Thermal profile and fluorescence reading	18
2	KIT CONTENT	4	11.3.2 Setup of samples/controls	19
3	STORAGE AND STABILITY OF REAGENTS	5	11.4 Preparation of the reaction mix	20
4	PRECAUTIONS FOR USE	6	11.5 Analysis of results	21
5	SAFETY RULES	8	11.5.1 Verify the run	21
5.1	General safety rules	8	11.5.2 Interpretation of results	23
5.2	Safety rules concerning the kit	8	11.6 TROUBLESHOOTING	25
6	MATERIAL REQUIRED, BUT NOT PROVIDED	9	12 DEVICE LIMITATIONS	27
6.1	Reagents	9	13 DEVICE PERFORMANCE	27
6.2	Instruments	9	13.1 Analytical specificity	27
6.3	Disposables	10	13.2 Analytical sensitivity: detection limit	27
7	INTRODUCTION	11	13.3 Analytical sensitivity: linear range	28
8	TEST PRINCIPLE	12	13.4 Reproducibility	28
9	PRODUCT DESCRIPTION	14	13.5 Diagnostic specificity	29
10	SAMPLE COLLECTION, MANIPULATION AND PRETREATMENT	15	13.6 Diagnostic sensitivity	29
10.1	Cytological samples	15	13.7 Accuracy	29
10.2	Histological samples	15	14 REFERENCES	30
10.3	Blood (PBMC)	16	15 RELATED PRODUCTS	31
11	PROTOCOL	17		
11.1	DNA extraction	17		
11.2	Internal control	17		

1 PRODUCT INFORMATION

1.1 Intended use

The REALQUALITY RS-MYCO P kit is an IVD for detection of the DNA of *Mycoplasma pneumoniae* (MYCO P). If used in combination with the REALQUALITY RQ-MYCO P STANDARD, code RQ-48-ST, it allows the quantification of the bacterial DNA present in the sample.

The test is based on Real-Time PCR on DNA extracted from human clinical samples.

This *in vitro* diagnostic test for detection and quantification of *Mycoplasma pneumoniae* is an auxiliary device for diagnosis and monitoring of *Mycoplasma pneumoniae* infections. It is recommended to use this kit as indicated in the instructions herein.

This manual refers to the following product:

REALQUALITY RS-MYCO P

Kit for detection and quantification of *Mycoplasma pneumoniae* (MYCO P) by Real-Time PCR.

[This product is in accordance with directive 98/79/EC \(Annex III\) on *in vitro* diagnostic medical devices \(CE marking\).](#)

Contains all reagents needed for Real-Time PCR.

Code	Product	PKG
RQ-S47-48	REALQUALITY RS-MYCO P	48 tests
RQ-S47-96	REALQUALITY RS-MYCO P	96 tests

2 KIT CONTENT

BOX RG

STORE AT -30 °C TO -20 °C

DESCRIPTION	LABEL	TUBE (T) or CAP COLOR	24 tests	48 tests	96 tests
Mastermix containing PCR reagents	MYCO P Real time mix	Violet	1 × 540 µL	2 × 540 µL	4 × 540 µL

BOX PC

STORE AT -30 °C TO -20 °C

DESCRIPTION	LABEL	TUBE (T) or CAP COLOR	24 tests	48 tests	96 tests
Positive control MYCO P (DNA fragment of the <i>Mycoplasma pneumoniae</i> genome)	PC MYCO P	Violet	1 × 50 µL	1 × 50 µL	2 × 50 µL
Positive control BG (DNA fragment of the β-globin gene)	PC BG	Blue	1 × 50 µL	1 × 50 µL	2 × 50 µL
Internal control (DNA fragment of the β-globin gene)	IC		2 × 125 µL	4 × 125 µL	8 × 125 µL

3 STORAGE AND STABILITY OF REAGENTS

Each component of the kit must be stored at the conditions indicated on the label of each box:

Box RG	Store at -30 °C to -20 °C
Box PC	Store at -30 °C to -20 °C

If stored at the recommended temperature all reagents are stable until the expiration date on the box.

Avoid degradation of the MYCO P Real time mix! The mix should NOT undergo more than two freeze/thaw cycles. If performing runs with low numbers of samples, it is recommended to aliquot the reagent beforehand. The MYCO P Real time mix contains fluorescent molecules and should be stored protected from direct light.

In order to avoid degradation of the positive controls and the internal control do NOT let them undergo more than three freeze/thaw cycles. If performing runs with low numbers of samples, it is recommended to aliquot the controls beforehand.

4 PRECAUTIONS FOR USE

- The kit must be used only as an IVD and be handled by qualified technicians who are trained in techniques of molecular biology applied to diagnostics.
- Before using the kit read the user manual carefully and completely.
- Keep the kit protected from heat.
- Please pay particular attention to the expiration date on the label of each box. Do not use any part of the kit past the expiration date.
- The reagents present in the kit must be considered an undividable unit. Do not use them separately or in combination with reagents from other kits or lots.
- The MYCO P Real time mix must be thawed at room temperature before use. Mix the solution by inverting the tube several times, then centrifuge briefly. Do NOT vortex!
- The positive controls and the internal control must be thawed at room temperature before use. Then centrifuge briefly.
- Work quickly, particularly if preparing the reactions at room temperature. If possible work on ice or on a cooling block.

In case of any doubt concerning storage conditions, box integrity or application of the method, please contact the technical support team at AB ANALITICA:

laboratorio@abanalitica.it

For nucleic acid amplification the user has to take the following precautions:

- Use filter-tips.
- In order to avoid contamination store biological samples, extracted DNA, amplification product and the internal and positive controls included in the kit separate from the MYCO P Real time mix.
- Set up pre- and post-PCR areas. Do not share instruments or consumables (pipettes, tips, tubes, etc.) between those areas.

- Change gloves frequently.
- Wash the bench surfaces with 5 % sodium hypochlorite.

5 SAFETY RULES

5.1 General safety rules

- Wear disposable gloves when handling reagents and clinical samples. Wash hands after the procedure.
- Do not pipet by mouth.
- No known diagnostic method can ensure the absence of infective agents. Therefore, consider all clinical samples to be potentially infectious and handle them accordingly.
- All devices that come into contact with clinical samples must be considered contaminated and disposed of as such. In case of accidental spilling of samples, clean up with 10 % sodium hypochlorite.

Material you use to clean must be disposed of in special containers for contaminated products.

- Clinical samples, contaminated materials and products must be decontaminated before disposal.

It is recommended to use one of the following decontamination methods:

- a) immerse for 30 minutes in a solution of 5 % sodium hypochlorite (1 volume of 5 % sodium hypochlorite solution on 10 volumes of contaminated fluid).
- b) autoclave at 121 °C for at least 2 hours (ATTENTION! Do not autoclave solutions containing sodium hypochlorite!!).

5.2 Safety rules concerning the kit

The risks for use of this kit are related to the single components.

Dangerous components: **none**.

The Material Safety Data Sheet (MSDS) of this device is available upon request.

6 MATERIAL REQUIRED, BUT NOT PROVIDED

6.1 Reagents

- Reagents for DNA extraction;
- DNase- and RNase-free sterile water;
- For quantitative analysis: REALQUALITY RQ-MYCO P STANDARD, code RQ-48-ST.

6.2 Instruments

- Laminar flow cabinets
Use while preparing the amplification mix in order to avoid contamination. It is recommended to use a different laminar flow cabinet when adding the extracted DNA and the positive controls / quantification standards.
- Micropipettes (range: 0.5 - 10 µL, 2 - 20 µL, 10 - 100 µL, 20 - 200 µL, 100 - 1000 µL)
- Microcentrifuge (max. 12,000 – 14,000 rpm)
- Plate centrifuge (optional)
- Real-Time PCR instrument

This kit has been validated on:

- Applied Biosystems 7500 Fast / 7500 Fast Dx Real-Time PCR System (ABI 7500 Fast / 7500 Fast Dx – *Applied Biosystems*)
- Applied Biosystems 7300 Real-Time PCR System (ABI 7300 – *Applied Biosystems*)
- Applied Biosystems StepOne / StepOnePlus™ Real-Time PCR System (ABI StepOne / StepOnePlus – *Applied Biosystems*)
- LightCycler® 480 Real-Time PCR System version II (LC 480 – *Roche*)
- LightCycler® 2.0 Real-Time PCR System (LC 2.0 – *Roche*)
- Dx Real-Time System (Bio-Rad Dx – *Bio-Rad*)
- CFX96 Real-Time PCR Detection System (Bio-Rad CFX96 – *Bio-Rad*)

The kit can be used on instruments that allow 25 µL of reaction volume and can read the fluorescence of the fluorophores FAM and JOE. The JOE fluorescence can also be read in the channels designated for CY3, HEX etc. Compatibility of the device with other commercially available

instruments has been asserted. For further information on instrument compatibility please contact the technical support team at AB ANALITICA.

6.3 Disposables

- Talc-free disposable gloves
- Disposable sterile filter-tips (range: 0.5 - 10 µL, 2 - 20 µL, 10 - 100 µL, 20 - 200 µL, 100 - 1000 µL)
- 96-well plates for Real-Time PCR with adhesive optical film, 0.1 - 0.2 mL tubes with optical caps or glass capillaries

7 INTRODUCTION

Mycoplasma pneumoniae is the most common etiologic agent for primary atypical pneumonia (PAP) occurring in children and younger age groups (below age of 30).

Mycoplasma pneumoniae are strictly aerobic bacteria that attaches to the ciliated respiratory epithelium by means of a special adhesin. Attaching to the cilia of epithelial cells it induces ciliostasis (block of the cilium) leading to subsequent destruction of the epithelium. This causes irritation of the respiratory tract and coughing.

The infection is asymptomatic in most cases but may cause severe pneumonia involving hematologic and neurologic manifestations in immunosuppressed patients. Infection of the respiratory tract by *Mycoplasma pneumoniae* is also causative for conditions like tracheobronchitis, pharyngitis and asthma. As *Mycoplasma pneumoniae* is capable of infecting other organs it is known to also cause pathological conditions of the central nervous system, liver, pancreas, blood, skin and joints.

The bacteria are spread through respiratory droplet transmission (droplets of saliva discharged by coughing). In persons with an active infection the pathogen can also be found in the sputum. Presence of the bacteria also in nose and throat indicates a diffuse infection.

Diagnosis of a *Mycoplasma pneumoniae* infection is traditionally difficult due to the use of serological methods and cell cultures which require time. This used to prevent prompt application of an effective therapy. Moreover, these bacteria often prove difficult to cultivate in the laboratory and for this reason frequently are not recognized as origin of the disease.

Detection of *Mycoplasma pneumoniae* using Real-Time PCR represents the most reliable and fastest method of diagnosis today. In addition, the Real-Time PCR allows exact quantification of the pathogen load.

8 TEST PRINCIPLE

The PCR method (Polymerase Chain Reaction) was the first method of DNA amplification described in literature (Saiki RK et al., 1985). It can be defined as an *in vitro* amplification reaction of a specific part of DNA (target sequence) by a thermostable DNA polymerase.

This technique was shown to be a valuable and versatile instrument of molecular biology: its application contributed to a more efficient study of new genes and their expression and has revolutionized for instance the fields of laboratory diagnostics and forensic medicine.

The Real-Time PCR represents an advancement of this basic research technology, providing the possibility to determine the number of amplified DNA molecules (amplicons) during the polymerase chain reaction (PCR).

In the system at hand monitoring the amplicons is based on primers/probes labeled with fluorescent molecules. These probes contain a reporter fluorophore and a molecule (quencher) that blocks the reporter's specific fluorescence. Fluorescent emission of the reporter is determined by its distance to the quencher. As long as a probe is not bound to a target sequence reporter and quencher are in close proximity and the reporter's fluorescence is blocked. Upon binding to a target sequence quencher and reporter become separated and the reporter can emit fluorescent light which in turn can be detected.

Typically, the main part of a Real-Time PCR run consists of 30 – 50 amplification cycles. A thermocycler equipped with a corresponding detector can record the fluorescence events at each cycle, thus monitoring the reaction in "real time".

The cycle at which the amplicon-related fluorescence becomes clearly distinguishable from the background is specific for each reaction and is correlated to the initial concentration of the target sequence. This cycle is called threshold cycle (Ct). The Ct value is used to determine the initial target concentration with the help of a standard curve. Such a standard curve is created amplifying solutions with known concentrations of the target sequence (Fig. 1).




Fig. 1: Creating a standard curve using standards with known concentrations.

The main advantage of Real-Time PCR compared to conventional techniques of amplification is the possibility to perform a semi-automated amplification. This means, extra steps necessary to visualize the amplification product can be avoided and the risk of contamination by post-PCR manipulation is reduced.

9 PRODUCT DESCRIPTION

The REALQUALITY RS-MYCO P kit, code RQ-S47, is an IVD for detection of *Mycoplasma pneumoniae* by amplification of a fragment of the gene coding for P1 cytoadhesin.

If used in combination with the product REALQUALITY RQ-MYCO P STANDARD, code RQ-48-ST, it allows the quantification of bacterial DNA molecules in the sample by means of a four-point standard curve (10^2 to 10^5 copies of bacterial DNA per reaction).

The positive controls supplied in this kit contain DNA fragments that correspond to the amplified gene region. As such, these controls are not harmful for the user.

The kit is designed to use an internal control that allows to detect inhibition of the PCR reaction, monitor the extraction process as well as identify false-negative samples. The internal control (the β -globin gene) is amplified in multiplex with the target pathogen. In cellular samples the endogenous β -globin gene is amplified. For acellular specimens the internal control is added as recombinant DNA containing the respective β -globin gene region.

The kit includes a ready-to-use mastermix that contains all reagents needed for the PCR as well as the components listed below.

- ROX™ is an inert colorant that exhibits stable fluorescent properties throughout all amplification cycles. On some Real-Time PCR instruments (Applied Biosystems, Stratagene etc.) it is used for normalization, in order to compensate differences between wells due to pipetting errors or limitations of the instrument.
- The dUTP/UNG system prevents contamination from previous amplification runs. The dUTPs are used to incorporate uracil residues into the amplification product during the amplification session. At the beginning of each new run the UNG enzyme degrades any single or double stranded DNA containing uracil. This way any amplification products from former sessions are eliminated.

10 SAMPLE COLLECTION, MANIPULATION AND PRETREATMENT

The detection of *Mycoplasma pneumoniae* can be performed on naso- and oropharyngeal swabs and samples of sputum, BAL, bronchoaspirate, expectoration, lung biopsy tissue and PBMC.

The device was tested with DNA extract of nasopharyngeal swabs, sputum, saliva, expectoration samples, PBMC and FFPE (formalin-fixed paraffin-embedded) tissue.

10.1 Cytological samples

Expectoration samples may be liquefied adding N-acetyl-L-cystein (NAC, 1.5 %) or dithiothreitol (DDT, 0.1 %).

Use sterile, disposable containers with a screw cap for sample collection.

Store the sample at +2 °C to +8 °C and extract the nucleic acids within 48 hours. If extraction is not feasible within 48 hours store the samples at -30 °C to -20 °C.

10.2 Histological samples

Histological samples include fresh, frozen or formalin fixed paraffin embedded tissue samples.

Tissue from biopsies can be used fresh if prepared only a few minutes after the tissue sample was taken. Alternatively, shock freeze the samples in liquid nitrogen and store the frozen samples at -80 °C.

For preparing fresh or frozen samples (up to 50 mg) for analysis, proceed without delay with mechanical disruption of the tissue (e.g. using a sterile scalpel) and subsequent digestion with Proteinase K.

For preparation of FFPE samples (formalin-fixed paraffin-embedded), remove the paraffin and then perform enzymatic digestion.

For fixation of tissue, it is recommended to use a sodium or potassium salt buffered, pH-neutral (pH 7) fixation solution with 10 % formalin (according to Lillie). Tissue that was fixed in non-buffered formalin or in Bouin's, Hollande's or acid based (i.e. Osmic acid) fixatives is not suited for DNA extraction. These solutions facilitate cross-link formation in the tissue, rendering it not digestible.

10.3 Blood (PBMC)

Sample collection should follow common routine, respecting all the usual sterility precautions (e.g. transport in sterile boxes without transport medium). The blood must be treated with EDTA. Other anticoagulation agents, like heparin, are strong inhibitors of the Taq polymerase and may impair the PCR. Store fresh blood at +2 °C to +8 °C and extract the nucleic acids within 4 hours. If extraction is not feasible within 4 hours store the samples at -30 °C to -20 °C.

Analysis is usually performed on a leukocyte pellet. Use density gradient separation (e.g. Ficoll-Hypaque).

11 PROTOCOL

11.1 DNA extraction

For DNA extraction AB ANALITICA recommends the QIAamp DNA Mini Kit (QIAGEN, Hilden, Germany). Refer to the manufacturer's manual for instructions and protocols for the different sample types.

This IVD device can be used with DNA extracted with the most common manual and automated extraction methods.

For further information regarding the compatibility of the device with different extraction methods please contact the technical support at AB ANALITICA.

11.2 Internal control

The kit includes an internal control consisting of a recombinant DNA fragment of the β -globin gene (BG). Use of this control is recommended for analysis of acellular samples. It allows to verify the extraction procedure and detect inhibition of the PCR.

The standardization of the internal control was performed adding 10 μ L of internal control to the sample and eluting in a volume of 60 μ L.

If the extraction system uses a different final elution volume, adjust the volume of internal control to be added to the sample, accordingly.

For correct use of the internal control follow the instructions provided by the extraction system manufacturer.

For any further information please contact the technical support at AB ANALITICA.

11.3 Instrument programming

11.3.1 Thermal profile and fluorescence reading

Set up the following thermal profile in your instrument:

	Step	Repeats	Time	(°C)
UNG Activation	1	1	02:00	50.0
Taq Activation	2	1	10:00	95.0
Amplification cycles	3	45	00:15 01:00	95.0 60.0 *

* Fluorescence detection step

The fluorophores to be read are:

- FAM for MYCO P
- JOE for BG

Select the two detection channels on your Real-Time PCR instrument:

	Name	Reporter Dye	Quencher Dye
ABI 7500 Fast / 7500 Fast Dx *	MYCO P	FAM	None
ABI 7300 *	BG (Internal Control)	JOE	None
ABI StepOne / StepOnePlus *			

	Name	Fluorophore	Filter
LC 480	MYCO P	FAM	465 - 510
	BG (Internal Control)	JOE	533 - 580

	Name	Fluorophore	Channel
LC 2.0	MYCO P	FAM	530
	BG (Internal Control)	JOE	560

	Name	Fluorophore
Bio-Rad Dx Bio-Rad CFX96	MYCO P	FAM
	BG (Internal Control)	JOE

* For instruments that require a passive reference (e.g. Applied Biosystems, Stratagene) make sure to select ROX for all wells in use.

Set the final reaction volume.

11.3.2 Setup of samples/controls

Set up samples, control(s) and standards (if needed) in the instrument software. Name each sample, control and standard accordingly.

Be careful to use the same position/order for your samples, controls and standards as for your real samples.

If you want to perform a quantitative analysis, enter the concentrations of the MYCO P standards (10^2 , 10^3 , 10^4 and 10^5 bacterial genome copies/reaction [c/rx]).

If you prefer to obtain quantification results in *bacterial genome copies per mL of clinical sample* (c/mL) you need to calculate and enter the standard concentrations depending on your extraction parameters. See the table below for examples.

STANDARD		Extraction Parameters		
		Example 1 $V_i = 200 \mu\text{L}$ $V_e = 50 \mu\text{L}$	Example 2 $V_i = 200 \mu\text{L}$ $V_e = 60 \mu\text{L}$	Example 3 $V_i = 200 \mu\text{L}$ $V_e = 100 \mu\text{L}$
STANDARD 1	10^2 c/rx	$5 \times 10^3 \text{ c/mL}$	$6 \times 10^3 \text{ c/mL}$	10^4 c/mL
STANDARD 2	10^3 c/rx	$5 \times 10^4 \text{ c/mL}$	$6 \times 10^4 \text{ c/mL}$	10^5 c/mL
STANDARD 3	10^4 c/rx	$5 \times 10^5 \text{ c/mL}$	$6 \times 10^5 \text{ c/mL}$	10^6 c/mL
STANDARD 4	10^5 c/rx	$5 \times 10^6 \text{ c/mL}$	$6 \times 10^6 \text{ c/mL}$	10^7 c/mL

V_i = initial volume (clinical sample volume used for extraction)

V_e = elution volume

11.4 Preparation of the reaction mix

Thaw the MYCO P Real time mix. After thawing homogenize the mix by inverting the tube several times. Do not vortex! Centrifuge briefly.

Work rapidly. If possible work on ice or a cooling block and in an area protected from direct light.

Note: A quantitative analysis requires the REALQUALITY RQ-MYCO P STANDARD, code RQ-48-ST.

Pipet 20 μL of the MYCO P Real time mix into the corresponding positions (wells of PCR plate, tubes, capillaries etc). Make sure to prepare sufficient positions for all samples, positive controls / quantification standards as well as the negative control.

Pipet 5 μL of extracted DNA, of negative control (sterile H_2O) or of positive controls / quantification standards into the corresponding positions.

Note: Thaw, mix and spin down the controls/standards before use.

Make sure no air bubbles remain in the wells/tubes/capillaries and centrifuge at 4000 rpm for 1 minute.

Load the samples on the instrument making sure to position/load the plate / tubes / capillaries correctly.

11.5 Analysis of results

After the PCR run is finished, view the analysis graph in logarithmic scale. Analyze the amplification results separately for *Mycoplasma pneumoniae* (MYCO P) and the β -globin gene (BG). Proceed as follows:

11.5.1 Verify the run

Before interpreting the results of the clinical samples you need to verify the PCR run. Evaluate the controls and/or standard curve according to the tables below.

A. QUALITATIVE ANALYSIS: Evaluating the controls

	RESULT	INTERPRETATION
Positive control BG	Amplification signal in JOE	Control and PCR worked correctly
	No amplification signal in JOE	No amplification of <i>BG</i> gene repeat the analysis
Positive control MYCO P	Amplification signal in FAM	Control and PCR worked correctly
	No amplification signal in FAM	No amplification of <i>Mycoplasma pneumoniae</i> DNA repeat the analysis
Negative control	Amplification signal in FAM and/or JOE	Contamination repeat the analysis
	No amplification signal in any channel	Control and PCR worked correctly

Only if all controls worked correctly, the run is suited for analysis.

B. QUANTITATIVE ANALYSIS: Evaluating controls and standard curve

	RESULT	INTERPRETATION
Positive control BG	Amplification signal in JOE	Control and PCR worked correctly
	No amplification signal in JOE	No amplification of <i>BG</i> gene repeat the analysis
Negative control	Amplification signal in FAM and/or JOE	Contamination repeat the analysis
	No amplification signal in any channel	Control and PCR worked correctly

INSTRUMENT	STANDARD CURVE PARAMETERS
ABI 7300 ABI 7500 Fast / 7500 Fast Dx ABI StepOne / StepOnePlus Bio-Rad Dx Bio-Rad CFX96	-3.60 < slope < -3.10 R ² > 0.99
LC 480	-3.60 < slope < -3.10
LC 2.0	1.8 < Efficiency < 2.1

The run is suited for analysis/interpretation, if all controls worked correctly and the standard curve parameters are in the specified range.

11.5.2 Interpretation of results

If the controls show the expected results, continue with the interpretation of the sample results. See the table below.

Target BG	Target MYCO P	INTERPRETATION
Ct (of amplification) < 34 (for whole blood: Ct < 32) §	Amplification signal	Positive for <i>Mycoplasma pneumoniae</i>
	No amplification signal	Negative for <i>Mycoplasma pneumoniae</i>
No amplification signal or Ct (of amplification) > 34 (for whole blood: Ct > 32) §	Amplification signal	* Positive for <i>Mycoplasma pneumoniae</i>
	No amplification signal	Not suitable for analysis Repeat DNA extraction

§ These Ct values refer to a DNA extraction from 200 µL of clinical sample and a final elution volume of 50 µL.

* **ATTENTION!** This assay has been optimized to favor amplification of the pathogen DNA. Therefore, the amplification signal of the control gene (β-globin, JOE fluorescence) may be delayed or absent in samples positive for *Mycoplasma pneumoniae*.

If quantification standards were included in the amplification run, the absolute number of *Mycoplasma pneumoniae* genome copies in the samples can be determined.

The exact number of viral genome copies can be determined only for results that are in the linear range of the device. See the table below for correct interpretation:

Quantification result <i>Mycoplasma pneumoniae</i>	INTERPRETATION (bacterial genome copies / reaction)
quantification result > 10 ⁷ copies/reaction	more than 10 ⁷
5 < quantification result < 10 ⁷ copies/reaction	exact quantity = quantification result
quantification result < 5 copies/reaction	less than 5

The pathogen load can be calculated as *genome copies/mL of clinical sample*, if the specific extraction parameters are included in the calculation. See the table below for examples.

Extraction parameters	Quantification result <i>Mycoplasma pneumoniae</i>	INTERPRETATION (bacterial genome copies / mL of clinical sample)
Example 1 V _i = 200 µL V _e = 50 µL	quantification result > 5 × 10 ⁸ copies/mL	more than 5 × 10 ⁸
	250 < quantification result < 5 × 10 ⁸ copies/mL	exact quantity = quantification result
	quantification result < 250 copies/mL	less than 250
Example 2 V _i = 200 µL V _e = 60 µL	quantification result > 6 × 10 ⁸ copies/mL	more than 6 × 10 ⁸
	300 < quantification result < 6 × 10 ⁸ copies/mL	exact quantity = quantification result
	quantification result < 300 copies/mL	less than 300
Example 3 V _i = 200 µL V _e = 100 µL	quantification result > 10 ⁹ copies/mL	more than 10 ⁹
	500 < quantification result < 10 ⁹ copies/mL	exact quantity = quantification result
	quantification result < 500 copies/mL	less than 500

V_i = initial volume (clinical sample volume used for extraction)
V_e = final elution volume

11.6 TROUBLESHOOTING

No amplification signals for positive controls / standards and samples

- *The instrument was not programmed correctly*
 - ▶ Repeat the amplification taking care of the instrument programming. Pay particular attention to the thermal profile, the fluorophores selected and that the positions of samples in the instrument setup correspond to the actual positions/order of the samples / controls / standards.
- *The reaction mix did not work correctly*
 - ▶ Make sure to store the MYCO P Real time mix at -30 °C to -20 °C. Avoid unnecessary freeze/thaw cycles. Store the mix protected from light.
 - ▶ Do not use the product past the expiration date reported on the label.

Very weak amplification signal for positive controls / standards

- *The positive controls / standard solutions were not stored correctly and have degraded*
 - ▶ Make sure to store the positive controls / standard solutions at -30 °C to -20 °C and do not let them undergo more than three freeze/thaw cycles.
 - ▶ Do not use the product past the expiration date.

Amplification signal of β -globin is very delayed or absent in the extracted sample (negative for *Mycoplasma pneumoniae*)

- *The extracted DNA was not suited for PCR and the reaction was inhibited*
 - ▶ Make sure to extract the nucleic acids correctly.
 - ▶ If an extraction method uses wash steps with solutions containing ethanol, make sure no ethanol residue remains in the DNA extract.
 - ▶ Use the extraction systems recommended in paragraph 11.1.
- *The clinical sample is not suited for analysis*
 - ▶ Make sure to correctly store and pretreat the clinical sample before performing the analysis.

In case of any further problems, please contact the technical support team at AB ANALITICA:

laboratorio@abanalitica.it

fax (+39) 049-8709510

tel. (+39) 049-761698

12 DEVICE LIMITATIONS

The kit can have reduced performance if:

- The clinical sample is not suitable for this analysis,
- The DNA is not suitable for PCR (due to the presence of PCR inhibitors or to the use of an inappropriate extraction method),
- The kit was not stored correctly.

13 DEVICE PERFORMANCE

The performance reported below is verified for all instruments this device has been validated on. For further information contact the technical support team at AB ANALITICA.

13.1 Analytical specificity

The specificity of the REALQUALITY RS-MYCO P kit, code RQ-S47, is guaranteed by an accurate and specific selection of primers and probe and by the use of stringent amplification conditions.

Alignment of primers and probe in the most important databases showed no non-specific pairing.

In order to analyze possible cross-reactions of this assay, samples positive for potentially cross-reactive pathogens were tested with this IVD. None of the tested pathogens gave a positive result.

13.2 Analytical sensitivity: detection limit

Serial dilutions of a quantification standard, ranging from 2 to 0.05 copies of bacterial genome/ μL , were tested in three consecutive experiments. Five microliters (5 μL) of each dilution were amplified in eight replicates per run and in multiplex with the internal control.

The results were analyzed using Probit analysis, as illustrated in Fig. 2.

The limit of the analytical sensitivity for the REALQUALITY RS-MYCO P kit ($p = 0.05$) on the ABI 7500 Fast Dx system is 1 bacterial genome copy/ μL DNA extract.




Fig. 2: Probit analysis for determination of the analytical sensitivity of the REALQUALITY RS-MYCO P kit (Applied Biosystems 7500 Fast Dx Real-Time PCR system). Displayed as bacterial genome copies / reaction.

13.3 Analytical sensitivity: linear range

The linear range of this assay was determined using a panel of dilutions of the quantification standards. Analysis was performed using linear regression. The linear range of the REALQUALITY RS-MYCO P kit on the ABI 7500 Fast Dx system is 5 to 10^7 bacterial genome copies / reaction.

13.4 Reproducibility

In order to determine the intra-assay variability (variability in one analysis session among replicates of the same sample) a dilution of 50 bacterial genome copies/ μL of the quantification standard (corresponding to a final amount of 250 copies/reaction) was amplified in eight replicates in one run.

The intra-assay variability coefficient of the method concerning the cycle threshold value (Ct) is 0.12 % on the ABI 7500 Fast Dx system.

In order to determine the inter-assay variability (variability over different analysis sessions of replicates of the same sample) the least concentrated quantification standard (20 bacterial genome copies/ μL) was amplified in duplicates in three consecutive runs. For each run, the variability coefficient was calculated from the Ct of the samples.

The inter-assay variability coefficient was calculated as the average of the variability coefficients for each run. The inter-assay variability coefficient on the ABI 7500 Fast Dx system is 0.33 %.

13.5 Diagnostic specificity

A statistically significant number of samples negative for *Mycoplasma pneumoniae* were tested simultaneously with the REALQUALITY RS-MYCO P kit and another CE IVD device or a reference method. From the obtained results the diagnostic specificity was calculated. The diagnostic specificity of this device is 100 %.

13.6 Diagnostic sensitivity

A statistically significant number of samples positive for *Mycoplasma pneumoniae* were tested simultaneously with the REALQUALITY RS-MYCO P kit and another CE IVD device or a reference method. From the obtained results the diagnostic sensitivity was calculated. The diagnostic sensitivity of this device is 100 %.

13.7 Accuracy

The accuracy was calculated as the ratio of the number of correct test results to the total number of executed tests. The accuracy of the REALQUALITY RS-MYCO P kit is 100 %.

14 REFERENCES

Bohte R et al. Thorax 50, 543–547, 1995

Jacobs E. Clin. Infect. Dis. 17(Suppl. 1), S79–S82, 1993

Jokinen C et al. Clin Infect Dis 32, 1141–1154, 2001

Saiki RK et al. Science 230, 1350-1354, 1985

Talkington DF et al. p. 57-84. in W. M. Scheld, W. A. Craig, and J. M. Hughes (ed.), Emerging Infections 5. American Society for Microbiology, Washington, D.C, 2001

Waites KB and Talkington DF. Clin Microbiol Rev. 17, 697–728, table of contents. doi: 10.1128/CMR.17.4.697-728, 2004

15 RELATED PRODUCTS

REALQUALITY RQ-MYCO P STANDARD

Ready-to-use quantification standards for quantification of *Mycoplasma pneumoniae* (MYCO P).

This product is in accordance with directive 98/79/EC (Annex III) on *in vitro* diagnostic medical devices (CE marking).

Code	Product	PKG
RQ-48-ST	REALQUALITY RQ-MYCO P STANDARD	4 × 60 µL



AB ANALITICA srl

Via Svizzera 16 - 35127 PADOVA, (ITALY)

Tel +39 049 761698 - Fax +39 049 8709510

e-mail: info@abanalitica.it